



## Sharps Container 2L Square with Flange and Seal

€8.81

In stock  
Only 999 left  
SKU 9970

 Add to Wish List

 Add to Compare

### Details

Sharps Container - 2L Square with Flange and Seal

Type - Inter-stacking (nesting), seamless, injection molded

Capacity to fill line - 1.7L

Capacity to spill line - 2L

Lid aperture - 90mm

Final Closure (red) - 90mm

Finger guard (yellow) - 90mm

Ideal application - Bedside, trolley clip-on or wall mounting; sharps, small waste eg band aids, needles, broken glass, swabs

Shape Square; with taper top to base

Dimensions - 190d x 130 w x 130 h mm



---

### **Important**

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other qualified healthcare provider. Consult your healthcare provider immediately if you think that you may have a medical emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified healthcare provider before starting any new treatment or with any questions that you may have regarding a medical condition or injury.

### **Related Products**

Check items to add to the cart or select all