



## Traverse Gazelle Basket Stretcher - Steel Mesh

€812.47

In stock  
Only **1000** left  
**SKU** 99166

 [Add to Wish List](#)

 [Add to Compare](#)

### Details

The Traverse Gazelle Basket Stretcher is a basic rescue stretcher designed to be used in a raise/lower operation. Its frame geometry, construction and materials provide its impressive strength. The top and lower rails are made using a 19mm robust steel tube. Durathene netting allows water and air to easily pass through the patient compartment.

Dimensions: 211 l x 62 w

Weight: 13kgs

Load rating: 408kgs

### Important

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other qualified healthcare provider. Consult your healthcare provider immediately if you think that you may have a medical



---

emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified healthcare provider before starting any new treatment or with any questions that you may have regarding a medical condition or injury.

**Related Products**

Check items to add to the cart or select all