





# Vitalograph Alpha 4 Spirometer

€3,433.73

In stock  
Only **1000** left  
**SKU** 99641

 Add to Wish List  
 Add to Compare



## Details

Vitalograph Alpha 4 Spirometer.

- Dedicated all-in-one desktop unit
- Simple to use
- Rapid and low cost testing.
- Clear and bright color screen.
- Accurate Fleisch pneumotachograph.
- Measures VC, FVC, FEV1, FEV1%, PEF, FEF 25-75 and more.
- On-screen visualisation of volume/time and flow/volume curves.
- Saves multiple pre-tests for subsequent post-testing.
- Compact Flash removable memory card slot.
- Presents the test information in clear, easy to read format.
- Links to PC, via USB, with report generation software.
- Fast, silent printer that prints 140 tests per roll
- Battery and mains power with integrated battery for off site testing
- Real time clock
- Autoclaveable flowhead
- Facility for user accuracy checks



- Automatically selects the best test
- Complete with stylish carry case

### **Important**

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other qualified healthcare provider. Consult your healthcare provider immediately if you think that you may have a medical emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified healthcare provider before starting any new treatment or with any questions that you may have regarding a medical condition or injury.

### **Related Products**

Check items to add to the cart or select all