



EgoPsoryl TA 30g

€10.10

In stock
Only **461** left
SKU 11515

 Add to Wish List

 Add to Compare

EgoPsoryl TA is a soothing corticosteroid-free gel, formulated to help control stable, moderate psoriasis and persistent dermatitis.

Details

EgoPsoryl TA is a soothing corticosteroid-free gel, formulated to help control stable, moderate psoriasis and persistent dermatitis.

Benefits

- Aids in the control of stable, moderate psoriasis (including psoriasis of the scalp) and persistent dermatitis.
- Soothing gel formulation calms the skin and relieves itching.
- Non staining.

Directions

Apply to affected skin 2 – 4 times daily. Avoid soap as soap may irritate skin. Cleanse with Pinetarsol, an anti-itch cleansing treatment which helps reduce inflammation and remove scales. For enhanced effect, apply petroleum jelly over EgoPsoryl TA.

Scalp: Apply at night if desired leave on during the day.

Cleanse with 'Sebitar' when needed.

Sensitive skin (eg face, children): Use Egoderm Cream for first 7 days. Apply EgoPsoryl TA over Egoderm Cream for 3 days. Then use EgoPsoryl TA alone.

Ingredients

Coal Tar Solution 5%, Phenol 0.5% and Sulfur 0.5%

Contains alcohol.