



QV Intensive Cleanser - 250g

€12.27

In stock
Only **500** left
SKU 10177

♥ Add to Wish List

▮ Add to Compare

A body wash formulated with 60% Glycerin to help maintain moisture whilst cleansing extremely dry and sensitive skin.



Details

QV Intensive Moisturising Cleanser is a rich, foaming body wash designed for extremely dry and sensitive skin. Scientifically formulated with 60% glycerin, the highest percentage in the QV range of washes and cleansers, the elevated percentage of moisturising ingredients allows for a reduction in water content and, because of this, the removal of preservatives.

Benefits:

- Free from preservatives, fragrance, colour and lanolin to suit sensitive skin types.
- Suitable for eczema, dry, sensitive and itchy skin.
- Helps to maintain skin's protective properties and doesn't strip away natural oils.
- Helps replace moisture to the skin as it cleanses.
- Vegan Friendly - contains no animal derived ingredients.



Directions:

For best results use this moisturising body wash daily in the shower. Massage cleanser onto your wet body, rinse and pat dry.

Cautions:

Avoid eyes. If cleanser enters eyes, rinse with clean water immediately. If irritation occurs discontinue use. Store below 30°C.

Ingredients:

Glycerin, Aqua (Water), Sodium Lauroyl Sarcosinate, Disodium Cocoamphodiacetate, Lauryl Betaine, Sodium Cocoyl Isethionate, Acrylates / C10-30 Alkyl Acrylate Crosspolymer, Styrene/Acrylates Copolymer, Aminomethyl Propanol.