



## QV Intensive Moisturiser - 450g

€20.94

In stock  
Only **500** left  
**SKU** 10183



♥ Add to Wish List

▮ Add to Compare

A highly emollient water-free soothing ointment for extremely dry, cracked sensitive skin. Helps to rehydrate and seal moisture in the the skin.



### Details

QV Intensive Moisturiser is a highly emollient balm that has been formulated to care for extremely dry or cracked skin. The moisturising balm helps to seal moisture into the skin, and the water-free formulation won't sting on cracked skin.

### Benefits:

- Water-free formulation won't sting when applied to inflamed or cracked skin.
- Suitable for dry, sensitive, itchy skin and eczema.
- Low-irritant formulation, free from fragrance, colour and lanolin to suit sensitive skin types.
- Helps to maintain skin's protective properties.
- Vegan Friendly - contains no animal derived ingredients.
- Body moisturising cream suitable for use with sensitive and eczema prone skin.



## Directions:

Apply body moisturiser liberally to the face and body and massage in gently. Use after bathing and showering.

## Cautions:

Avoid eyes. For external use only. Store below 30°C.

## Ingredients:

Light Liquid Paraffin (Paraffinum Liquidum), Petrolatum, Isopropyl Myristate, Polyethylene, Cetearyl Alcohol, Silica, Dimethicone