



AEROSWAB 15cm Single End Applicators Bag/100

€3.79

In stock
SKU ACA150


 Add to Wish List
 Add to Compare
Single End Applicator



Details

Size: 15cm **Pricing Unit:** bag/100 **Carton QTY:** 240 The AEROSWAB is a high quality cotton swab and is ideal for cleansing wounds when used in conjunction with an antiseptic or saline solution.

AEROSWAB Cotton Tips

AEROSWAB cotton tip applicators are manufactured from premium quality cotton and attached to sturdy applicators. 



More Information

| | |
|------------------------|----------|
| MANUFACTURER_MS | Germany |
| MODEL NUMBER_MS | K61SD521 |
| ECOLABEL_MS | yes |