



## AEROCREPE Elastic Crepe Bandage 10cm x 4M Wrap/12

€38.02

In stock  
SKU ACE100

 Add to Wish List

 Add to Compare

Elastic Crepe Bandage 10cm




### Details

**Size:** 10cm x 4M **Pricing Unit:** wrap/12 **Carton QTY:** 20 The AEROCREPE Bandage is a high quality cotton crepe bandage, recommended for low support requirements and can be used for the retention of dressings.

### Elastic Crepe

AEROCREPE elastic cotton crepe bandages are of an elastic conforming nature which is durable and suitable for immobilisation and compression techniques.

### Easy to Open Wrapper

AEROCREPE elastic cotton crepe bandages have a easy open tab built into the wrapper of the bandage. 



### More Information

|                        |          |
|------------------------|----------|
| <b>MANUFACTURER_MS</b> | Germany  |
| <b>MODEL NUMBER_MS</b> | K61SD521 |
| <b>ECOLABEL_MS</b>     | yes      |