



## LED Magnifying Lamp with Table Clamp (12cm diameter, 115cm extension)

€1,083.47

In stock  
SKU ACR05-ML

 Add to Wish List


 Add to Compare

Magnifying Lamp with Table Clamp

### Details

**Size:** 12cm diameter **Pricing Unit:** unit/1 **Carton QTY:** 2

### Magnifying Lamp

A magnifying lamp is integral to accurately diagnose and treat small injuries and provide sufficient light for casualty treatment. The AeroSupplies LED Magnifying Lamp comes with a table clamp and has a 1200mm reach/extension. It has lens 6 diopter magnification (2.5x) and is has 10 watt LED lights (6000K). You can purchase a mobile base ([ACR05-MB](#)) or wall mount ([ACR05-WM](#)) separately that are both compatible with this magnifying lamp for different applications. 



### More Information

<b>MANUFACTURER_MS</b>	Germany
<b>MODEL NUMBER_MS</b>	K61SD521
<b>ECOLABEL_MS</b>	yes