




## Small Stainless Steel Trolley with Drawer 50 x 50 x 97cm

€2,560.35

In stock  
**SKU** ACR12

 Add to Wish List

 Add to Compare  
Trolley with Drawer



### Details

**Size:** 50 x 50 x 97cm **Pricing Unit:** unit/1 **Carton QTY:** 1

### Medical/Instrument Trolleys

These units provide a sturdy surface and storage for medical supplies of instruments. Stainless steel construction makes for a durable, aesthetic and easy to disinfect finish. Supplied flat packed with minor assembly required.

### Product Features:

- 1 Drawer: 125mm deep
- High quality 304 grade stainless steel
- Safety fold to all edges
- Rubber wheels, front 2 locking
- Full extension telescopic ball bearing drawer runners



### More Information

<b>MANUFACTURER_MS</b>	Germany
<b>MODEL NUMBER_MS</b>	K61SD521
<b>ECOLABEL_MS</b>	yes