



AEROCLEANSE Sanitiser Station suitable for 1L bottle

€1,046.71

In stock
SKU ACSS01

 Add to Wish List

 Add to Compare

AeroCleanse Antibacterial Hand Gel

Details

Zero Touch, Zero Contamination

The patented AeroCleanse™ Sanitising Station is designed to prevent cross-contamination. This no-touch, foot-controlled, dispensing station is constructed from durable powder-coated aluminum with a steel base. The station is built-to-last and dispensing properly, every time guaranteed.

Technical Specifications

- Durable powder-coated steel
- 100% Mechanical
- No Contact with Hands
- No Batteries
- Height: 100cm





More Information

MANUFACTURER_MS	Germany
MODEL NUMBER_MS	K61SD521
ECOLABEL_MS	yes