



AEROPAD Non-Adherent Eye Pads 5.5cm x 7.7cm Box/50

€31.57

In stock
SKU AEP1S

 Add to Wish List

 Add to Compare

Low Adherent Eye Pad



Details

Size: 5.5 x 7.7cm **Pricing Unit:** box/50 **Carton QTY:** 24 The AEROPAD Non-Adherent dressing is recommended for treating cuts, grazes and minor burns using specialized SupaSorb technology.

Premium Absorbency

- Premium absorbent inner material for weeping and bleeding wounds

Non-adherent

- Special Non-adherent film on each side of the pad minimizing wound trauma when changing dressings.

Individually packed and sterile

- Dressings are available individually on demand and sterile for wound cleanliness.





More Information

MANUFACTURER_MS	Germany
MODEL NUMBER_MS	K61SD521
ECOLABEL_MS	yes