



RESPONDER 4 Series Plastic Rugged First Aid Kit 23.3 x 26.6 x 9.8cm

€193.56


In stock
SKU AFAK4C

 Add to Wish List

 Add to Compare

RESPONDER 4 Series Plastic Rugged First Aid Kit 23.3 x 26.6 x 9.8cm

Details

Ensure workplace compliance for first aid with this comprehensive range, containing a wide scope of products to cater for any workplace emergencies. Respond promptly to incidents with a variety of outer cases to suit all types of environments. Internal compartments and shelves keep product organised, and all plastic and metal case options are able to be wall-mounted. The soft pack style is ideal for mobile workers, while the durable plastic cases are recommended for industrial areas. **Contents meet Australian Standard AS2675-1983** This kit is available in 5 case options: a Robust rugged plastic case, Versatile Softpack reinforced fabric bag, a Weather-resistant plastic case and a Neat plastic case all our cases and bags are available with custom labelling. 

This First Aid Kit comes fully stocked and as part of our **Instant Restock™** system, you can refill and service your client's kits quickly with refill pouches designed to fit into each case type. They provide quick access to all product through a unique centre slit pocket without the need of unpacking each individual product.

KIT CONTENTS:

CODE	DESCRIPTION	QUANTITY
AP10125-I	AEROPLAST Plastic Plaster 72x19mm	25 pc
AFP50125-I	AEROPLAST Fabric Premium Strip 72mm x 19mm	25 pc
AB3-I	AEROBURN Burn Gel Sachet 3.5g	5 pc
AAS50-I	AEROAID Antiseptic Spray 50ml	1 pc



AFS001-I	AEROSHIELD CPR Face Shield	1 pc
AGNPF02-I	AEROGLOVE Nitrile Examination Gloves	5 Pairs
AGS753S-I	AEROSWAB Gauze Swabs 7.5cm x 7.5cm 3's	5 Pks
AW1000-I	AEROWASH Eye Wash Ampoule 15ml	8 pc
AW8001-I	AEROWIPE Cleansing Wipe	10 pc
ASP001-I	AEROPROBE Splinter Probes 3.7cm	10 pc
APD50S-I	AEROPAD Low Adherent Dressing 5cm x 5cm	6 pc
APD100S-I	AEROPAD Low Adherent Dressing 7.5cm x 10cm	3 pc
APD101S-I	AEROPAD Low Adherent Dressing 10cm x 10cm	2 pc
AF50-I	AEROFORM Conforming Bandage 5cm x 4M	3 Roll
AF75-I	AEROFORM Conforming Bandage 7.5cm x 4M	3 Roll
AFHP100-I	AEROFORM Heavyweight Conforming Bandage 10cm x 4m	1 Roll
AP325-I	AEROPORE Microporous Tape 2.5cm x 5M	1 Roll
AWD14S-I	AEROWOUND BPC Wound Dressing #14	1 Roll
AWD15S-I	AEROWOUND BPC Wound Dressing #15	1 Roll
ACD1020S-I	AEROWOUND Combine Dressing 10cm x 20cm	1 pc
ABC10-I	AEROBAND Triangular Bandage 110 x 110 x 155cm	2 pc
ATB130-I	AERORESCUE Emergency Rescue Blanket	1 pc
AEP1S-I	AEROPAD Eye Pads	4 pc
AII2500-I	AEROCOOL Instant Ice Pack 80g	1 pc
ASC13-I	AEROINSTRUMENTS Scissor 13cm	1 pc
AF13-I	AEROINSTRUMENTS Tweezer 11cm	1 pc
AAB100-I	AEROSUPPLIES Amputated Parts Bags Assorted Sizes	3 pc
ASN096P-I	AEROSUPPLIES Notebook & Pen	1 pc
ASP12-I	AEROPINS Safety Pin	12 pc
AGL001-I	AEROGUIDE First Aid Leaflet	1 pc



More Information

MANUFACTURER_MS	Germany
MODEL NUMBER_MS	K61SD521
ECOLABEL_MS	yes