



RESPONDER 4 Series Plastic Neat First Aid Kit Small 27 x 28 x 9.5cm

€196.76

In stock
SKU AFAK4P

 Add to Wish List

 Add to Compare


RESPONDER 4 Series Plastic Neat First Aid Kit Small 27 x 28 x 9.5cm

Details

Case: Plastic

Size: 27 x 28 x 10cm

Pricing Unit: unit/1

Ensure workplace compliance for first aid with this comprehensive range, containing a wide scope of products to cater for any workplace emergencies. Respond promptly to incidents with a variety of outer cases to suit all types of environments. Internal compartments and shelves keep product organised, and all plastic and metal case options are able to be wall mounted. The softpack style is ideal for mobile workers, while the durable plastic cases are recommended for industrial areas. **Contents meet Australian Standard AS2675-1983** This kit is available in 5 case options: a Robust rugged plastic case, Versatile Softpack reinforced fabric bag, a Weather resistant plastic case and a Neat plastic case - all our cases and bags are available with custom labelling. 

This First Aid Kit comes fully stocked and as part of our **Instant Restock™** system you can refill and service your client's kits quickly with refill pouches designed to fit into each case type. They provide quick access to all product through a unique centre slit pocket without the need of unpacking each individual product.

KIT CONTENTS:

| CODE | DESCRIPTION | QUANTITY |
|-----------|-----------------------------------|----------|
| AP10125-I | AEROPLAST Plastic Plaster 72x19mm | 25 pc |



| | | |
|------------|---|---------|
| AFP50125-I | AEROPLAST Fabric Premium Strip 72mm x 19mm | 25 pc |
| AB3-I | AEROBURN Burn Gel Sachet 3.5g | 5 pc |
| AAS50-I | AEROAID Antiseptic Spray 50ml | 1 pc |
| AFS001-I | AEROSHIELD CPR Face Shield | 1 pc |
| AGNPF02-I | AEROGLOVE Nitrile Examination Gloves | 5 Pairs |
| AGS753S-I | AEROSWAB Gauze Swabs 7.5cm x 7.5cm 3's | 5 Pks |
| AW1000-I | AEROWASH Eye Wash Ampoule 15ml | 8 pc |
| AW8001-I | AEROWIPE Cleansing Wipe | 10 pc |
| ASP001-I | AEROPROBE Splinter Probes 3.7cm | 10 pc |
| APD50S-I | AEROPAD Low Adherent Dressing 5cm x 5cm | 6 pc |
| APD100S-I | AEROPAD Low Adherent Dressing 7.5cm x 10cm | 3 pc |
| APD101S-I | AEROPAD Low Adherent Dressing 10cm x 10cm | 2 pc |
| AF50-I | AEROFORM Conforming Bandage 5cm x 4M | 3 Roll |
| AF75-I | AEROFORM Conforming Bandage 7.5cm x 4M | 3 Roll |
| AFHP100-I | AEROFORM Heavyweight Conforming Bandage 10cm x 4m | 1 Roll |
| AP325-I | AEROPORE Microporous Tape 2.5cm x 5M | 1 Roll |
| AWD14S-I | AEROWOUND BPC Wound Dressing #14 | 1 Roll |
| AWD15S-I | AEROWOUND BPC Wound Dressing #15 | 1 Roll |
| ACD1020S-I | AEROWOUND Combine Dressing 10cm x 20cm | 1 pc |
| ABC10-I | AEROBAND Triangular Bandage 110 x 110 x 155cm | 2 pc |
| ATB130-I | AERORESCUE Emergency Rescue Blanket | 1 pc |
| AEP1S-I | AEROPAD Eye Pads | 4 pc |
| AII2500-I | AEROCOOL Instant Ice Pack 80g | 1 pc |
| ASC13-I | AEROINSTRUMENTS Scissor 13cm | 1 pc |
| AF13-I | AEROINSTRUMENTS Tweezer 11cm | 1 pc |
| AAB100-I | AEROSUPPLIES Amputated Parts Bags Assorted Sizes | 3 pc |
| ASN096P-I | AEROSUPPLIES Notebook & Pen | 1 pc |
| ASP12-I | AEROPINS Safety Pin | 12 pc |
| AGL001-I | AEROGUIDE First Aid Leaflet | 1 pc |



More Information

| | |
|------------------------|----------|
| MANUFACTURER_MS | Germany |
| MODEL NUMBER_MS | K61SD521 |
| ECOLABEL_MS | yes |