



## OPERATOR 5 Series Metal Tough First Aid Kit 42 x 28.5 x 13cm

€374.16

In stock  
SKU AFAK5M

 Add to Wish List

 Add to Compare


OPERATOR 5 Series Metal Tough First Aid Kit 42 x 28.5 x 13cm

### Details

**Case:** Metal

**Size:** 42 x 29 x 13cm

**Pricing Unit:** unit/1

Compliance with workplace first aid regulations is achieved across all Australian states and territories with this broad kit range. Workplace injuries in higher risk areas can be treated effectively with the quality selection of First Aid supplies, complemented by a variety of outer cases to suit any scenario. Our metal cabinets can be wall mounted and key-lockable for security. For demanding environments, equip yourself with one of the robust plastic cases which are strong yet lightweight, optimising their portability. **Contents meet Australian Standard AS2675-1983**  This kit is available in 5 case options: a Robust rugged plastic case, Versatile Softpack reinforced fabric bag, a Weather resistant plastic case and a Neat plastic case - all our cases and bags are available with custom labelling.

This First Aid Kit comes fully stocked and as part of our **Instant Restock™** system you can refill and service your client's kits quickly with refill pouches designed to fit into each case type. They provide quick access to all product through a unique centre slit pocket without the need of unpacking each individual product.

### KIT CONTENTS:

CODE	DESCRIPTION	QUANTITY
AP10125-I	AEROPLAST Plastic Plaster 72x19mm	25 pc
AFP50125-I	AEROPLAST Fabric Premium Strip 72mm x 19mm	25 pc



AB3-I	AEROBURN Burn Gel Sachet 3.5g	5 pc
PTD06-I	AERODIAGNOSTIC Diagnostic Penlight	1 pc
AAS50-I	AEROAID Antiseptic Spray 50ml	1 pc
AFS001-I	AEROSHIELD CPR Face Shield	1 pc
AGNPF02-I	AEROGLOVE Nitrile Examination Gloves	5 Pairs
AGS753S-I	AEROSWAB Gauze Swabs 7.5cm x 7.5cm 3's	9 Pks
AW1000-I	AEROWASH Eye Wash Ampoule 15ml	8 pc
AW8001-I	AEROWIPE Cleansing Wipe	10 pc
ASP001-I	AEROPROBE Splinter Probes 3.7cm	10 pc
APD50S-I	AEROPAD Low Adherent Dressing 5cm x 5cm	6 pc
APD100S-I	AEROPAD Low Adherent Dressing 7.5cm x 10cm	3 pc
APD101S-I	AEROPAD Low Adherent Dressing 10cm x 10cm	2 pc
AF50-I	AEROFORM Conforming Bandage 5cm x 4M	3 Roll
AF75-I	AEROFORM Conforming Bandage 7.5cm x 4M	3 Roll
AFHP100-I	AEROFORM Heavyweight Conforming Bandage 10cm x 4.5m	2 Roll
AP325-I	AEROPORE Microporous Tape 2.5cm x 5M	1 Roll
AWD13S-I	AEROWOUND BPC Wound Dressing #13	2 Roll
AWD14S-I	AEROWOUND BPC Wound Dressing #14	3 Roll
AWD15S-I	AEROWOUND BPC Wound Dressing #15	3 Roll
ACD1020S-I	AEROWOUND Combine Dressing 10cm x 20cm	1 pc
ABC10-I	AEROBAND Triangular Bandage 110 x 110 x 155cm	4 pc
ATB130-I	AERORESCUE Emergency Rescue Blanket 127 x 180cm	1 pc
AEP1S-I	AEROPAD Eye Pads 5.5cm x 7.7cm	4 pc
AII2500-I	AEROCOOL Instant Ice Pack 80g	1 pc
ND90-I	AEROHAZARD Needle Disposal Unit 90ml Yellow	1 pc
TRYLD1-I	AEROHAZARD Disposable Dressing Tray	1 pc
PPC60-I	AEROHAZARD Portion Cup Plastic 60ml	1 pc
ASC13-I	AEROINSTRUMENTS Scissor 13cm	1 pc
AF13-I	AEROINSTRUMENTS Tweezer 11cm	1 pc
AAB100-I	AEROSUPPLIES Amputated Parts Bags Assorted Sizes	3 pc
ASN096P-I	AEROSUPPLIES Notebook & Pen	1 pc
ASP12-I	AEROPINS Safety Pins	12 pc
AGL001-I	AEROGUIDE First Aid Leaflet	1 pc



## More Information

<b>MANUFACTURER_MS</b>	Germany
<b>MODEL NUMBER_MS</b>	K61SD521
<b>ECOLABEL_MS</b>	yes