



AEROFIX Underwrap Tape 5cm x 10M

€7.94

In stock
Only **671** left
SKU AFR50



 Add to Wish List

 Add to Compare

Aerofix Non-Woven Fixation Tape



Details

Size: 5cm x 10M **Pricing Unit:** unit/1 **Carton QTY:** 100 AeroFix non-woven fixation tape is a highly conformable, latex-free solution for dressing retention and underwrap.

SupaBond adhesive

- SupaBond adhesive, exclusively engineered for unrivalled adhesion, provides a secure adhesion across the site, preventing lift up or edge curling.

Soft and Conformable

- A non-woven backing is smooth and soft, and adheres to body contours using SupaBond adhesion.

Contraindications / Precautions

- AeroFix should not be used on infected wounds unless under the supervision of a healthcare practitioner.
- As with all adhesive products, apply and remove carefully from sensitive or fragile skin.
- Seek medical advice immediately for all serious wounds and burns. If redness or discomfort occurs, discontinue use and consult a physician.

This product is not a substitute for direct medical attention from a healthcare practitioner. 



More Information

MANUFACTURER_MS	Germany
MODEL NUMBER_MS	K61SD521
ECOLABEL_MS	yes