



# AEROFORM Size A Newborn Elastic Tubular Bandage 4cm x 1M

€3.56

In stock  
SKU AGA01

 Add to Wish List

 Add to Compare

AEROFORM Size A Newborn Elastic Tubular Bandage 4cm x 1M


## Details

**Size:** 4cm x 1M **Pricing Unit:** unit/1 **Carton QTY:** 10 AEROFORM Tubular Bandages are an elasticised support bandage suited to treatment of sports injuries, dislocations and sprains. The Size A tubular bandage is suitable for use on infant limbs.

## Self Supporting

AEROFORM Tubular Bandages are held in place without tying or tapes, resulting in a comfortable long wearing bandage.

## Provides Even Support

The AEROFORM Tubular Bandages provides even support across damaged joints and muscles due to its elastic self supporting properties resulting in a comfortable support bandage. Available in 1m lengths or 10m rolls. 



### More Information

<b>MANUFACTURER_MS</b>	Germany
<b>MODEL NUMBER_MS</b>	K61SD521
<b>ECOLABEL_MS</b>	yes