





AEROSWAB Non-Sterile White Gauze Swab 10 x 10cm Pack/100 (8PLY)

€5.77

In stock
SKU AGS100

 Add to Wish List
 Add to Compare
8 ply Gauze Swab



Details

Size: 10 x 10cm **Pricing Unit:** pack/100 **Carton QTY:** 80 The AEROSWAB™ Gauze Swab is highly absorbent cotton gauze swab that when soaked in fluid becomes an ideal item to clean and prepare wounds for dressing.

Cotton Gauze

- AEROSWAB™ cotton gauze has a balanced combination of softness, abrasive texture and absorbency for effectively cleaning wounds.

Multiple Sizes Available

- AEROSWAB™ gauze swabs are available in a number of different sizes to suit any application.





More Information

| | |
|------------------------|----------|
| MANUFACTURER_MS | Germany |
| MODEL NUMBER_MS | K61SD521 |
| ECOLABEL_MS | yes |