





AEROSWAB Sterile White Gauze Swab 7.5 x 7.5cm Box/25 (Packs of 5)

€7.47

In stock
SKU AGS755S

-  Add to Wish List
 -  Add to Compare
- Gauze Swab - 5 Pcs



Details

Size: 7.5 x 7.5cm **Pricing Unit:** box/25 **Carton QTY:** 24 The AEROSWAB™ Gauze Swab is highly absorbent cotton gauze swab that when soaked in fluid becomes an ideal item to clean and prepare wounds for dressing.

Cotton Gauze

- AEROSWAB™ cotton gauze has a balanced combination of softness, abrasive texture and absorbency for effectively cleaning wounds.

Multiple Sizes Available

- AEROSWAB™ gauze swabs are available in a number of different sizes to suit any application.

Sterile

- AEROSWAB™ gauze swabs are also available in a sterile form for minimising infection risk when dressing clean wounds





More Information

MANUFACTURER_MS	Germany
MODEL NUMBER_MS	K61SD521
ECOLABEL_MS	yes