



AEROPAD Non-Adherent Dressing 7.5 x 10cm Box/50

€12.40

In stock
SKU APD100S

 Add to Wish List

 Add to Compare

Low Adherent Dressing Pad Box/50



Details

Size: 7.5 x 10cm **Pricing Unit:** box/50 **Carton QTY:** 16 The AEROPAD Non-Adherent dressing is recommended for treating cuts, grazes and minor burns using specialized SupaSorb technology.

Premium Absorbency

- Premium absorbent inner material for weeping and bleeding wounds

Non-Adherent

- Special Non-adherent film on each side of the pad minimizing wound trauma when changing dressings.

Individually packed and sterile

- Dressings are available individually on demand and sterile for wound cleanliness.





More Information

| | |
|------------------------|----------|
| MANUFACTURER_MS | Germany |
| MODEL NUMBER_MS | K61SD521 |
| ECOLABEL_MS | yes |