



AEROSHIELD Medium Nitrile Finger Cots Box/100

€10.83

In stock
SKU ASN100-M

♥ Add to Wish List

▬ Add to Compare
Nitrile Finger Cots



Details

Size: Medium **Pricing Unit:** box/100 **Carton QTY:** 100 AeroShield Nitrile Finger Cots are a high quality range of finger cots/stalls available in a range of sizes for first aid use.

Nitrile Finger Cots

Features

- Ambidextrous and anatomically shaped for comfort
- Powder free and optimum wet and dry grip
- Excellent flexibility and puncture strength
- Blue colour for easy detection of splits
- Pre-rolled for convenience
- Anti-static





More Information

MANUFACTURER_MS	Germany
MODEL NUMBER_MS	K61SD521
ECOLABEL_MS	yes