





## AEROWASH Drop Eyewash Bottle 100mL

€8.33

In stock  
**SKU** AWDREF100

-  Add to Wish List
  -  Add to Compare
- Drop Eyewash Refill



### Details

**Size:** 100ml **Pricing Unit:** unit/1 **Carton QTY:** 25 The innovative DROP eyewash station refill bottles all fit into the DROP eyewash station ([AWDS01](#)), which is a leader in innovation with a drop down dust cover, integrated mirror, integrated eye cup and a secure safety seal. DROP eyewash refills are a popular solution for convenient eye and wound irrigation.

### Strictly Sodium Chloride

- DROP Eyewash is made up strictly of Sodium Chloride 0.9% sterile solution which is safe and cleansing for eyes and wounds

### Secure Twist top opening

- DROP Eyewash refills have a secure twist top opening that keeps the saline sterile and secure but makes it easily available when necessary.

### Instructions

- Instructions for use are clearly indicated on the DROP Eyewash saline bottles

Also Available in [250ml](#) and [500ml](#) Refills



**More Information**

<b>MANUFACTURER_MS</b>	Germany
<b>MODEL NUMBER_MS</b>	K61SD521
<b>ECOLABEL_MS</b>	yes