



Polar H7 Bluetooth Smart Heart Rate Sensor



Out of stock
[Only 352 left](#)

Product availability

Product Name	Qty
Polar H7 Bluetooth Smart Heart Rate Sensor - black-S	88.0000
Polar H7 Bluetooth Smart Heart Rate Sensor - black-M	88.0000
Polar H7 Bluetooth Smart Heart Rate Sensor - blue-S	88.0000
Polar H7 Bluetooth Smart Heart Rate Sensor - blue-M	88.0000

SKU
 Polar_H7_Blueto
 oth

 Add to Wish List

 Add to Compare

Apart from the products which are on offered on the website, we have a well-trained staff who is ready to help you and to provide you with professional health care. See the full list of the hospitals where you can ask for a nurse and find the one which is near you.



Details

When it comes to caregiving, it's about creating comfort for those who are near us. We know it very well, because every single team member of our company works as a caregiver or used to work as one

Many of you face caregiving for the first time and find yourselves unprepared for this mission. Understanding this, we've created a regularly updated blog which guides you through the basics of caregiving and gives you the tips. Regardless of the impairment of the individual you must take care of, here you'll find the necessary information and assets which you need. Caregiving becomes much easier if you have appropriate equipment for that. Speaking of equipment - anything, literally anything that you will or may need to be a carer, you'll find here, in this online store.

If you are a fan of a specific brand, you can filter items using this criterion. If you are ready to consider all possible offers, please check out all categories. Everything for bath safety like bariatric bathroom aids, bathlifts, bathtub transfer benches, bath mats, cast protectors, commode chairs, shower chairs and bath benches, toilet safety frames and rails and so on. Everything related to a bedroom and daily living is also here. Dozens of models of home blood pressure monitors, electronic alerts systems, phones for seniors are available; dental supplies, chair cushions, bibs for adults, personal care items are at your disposal in this store. Apart from that, you'll find a huge number of equipment in Lift Chairs section, Medical, Nutrition, Walking Aids, Wheelchairs and even Scooters department.

If you or someone you love needs care, you may find a caregiver or a doctor near you. Just use our website to find a location and a number to find a doctor or a caregiver near you. The best stuff works in our store so that you can give and get perfect care.



More Information

BRAND_MS	Stander EZ
MANUFACTURER_MS	Italy
MODEL NUMBER_MS	K61SD521
TYPE_MS	Portable
SPECIFIC USE_MS	Hospital Bed
SHIPPING WEIGHT_MS	2,5 lb
OPERATING TEMPERATURE RANGE_MS	10°-30°C
LARGE SIZE CUFF	8.75"-16.5"
TOUCH SENSITIVE BUTTONS	yes
POWER_MS	(4) AA batteries
CASE MATERIAL	Metal
FEATURES_MS	Systolic blood pressure; Diastolic blood pressure; Pulse rate; Irregular Heartbeat; Detection and Hypertension Risk Indicator
ECOLABEL_MS	yes
WORKING HOURS_MS	30 minutes

Related Products

We found other products you might like!