




LIFEPAK CR2 Semi-Automatic Defibrillator with Wi-Fi

€2,596.67

In stock
SKU CR-2

 Add to Wish List


 Add to Compare

Semi-Automatic Defibrillator AED WiFi

Details

Size: - **Pricing Unit:** unit/1 **Carton QTY:** 1

- 8-year warranty
- Wireless connectivity
- Fastest time to 1st shock
- CPR Insight
- Durable
- 4 years standby battery life

Sudden cardiac arrest (SCA) can happen to anyone, anywhere. Immediate treatment is vital. A victim's chance of survival dramatically decreases for every minute without treatment. That's why public access defibrillators are so important: they put lifesaving technology where it can do the most good. So when an emergency happens, you should have nothing less than the best.  **Layered design** Layered design with easy to follow bold graphics. Both trained and untrained AED users clearly know how to begin. **QUIK-STEP™ electrodes** Peel directly off the base for faster placement. **cprINSIGHT analysis technology** Analyses for shockable rhythm during chest compressions with no need to pause. **Metronome and CPR coaching** Sets an effective pace and audibly guides users. **ClearVoice™ technology** Detects background noise and adjusts volume accordingly for clear instruction. **Fully automatic** Available in fully or semi-automatic models. **Highest available escalating energy** Up to 360J for more effective shocks as needed. **Bilingual** Toggle between two pre-set languages when using the device. **Child mode** Child mode delivers lower energy levels appropriate for young children without



having to change electrodes. **LIFEPAK TOUGH™** IP55 rating for challenging environments. **8-year warranty** Backed by an 8-year warranty. **LIFELINKcentral AED Program Manager** Monitor AED programs by tracking AED status, sending patient data to emergency responders and hospitals, detecting AED locations and other tools.



More Information

MANUFACTURER_MS	Germany
MODEL NUMBER_MS	K61SD521
ECOLABEL_MS	yes