




## CARDIACT Large Waterproof Tough AED Case 43 x 38 x 15.4cm

€245.55

In stock  
SKU CW-20


 Add to Wish List

 Add to Compare

Large Waterproof Tough AED Case



### Details

**Size:** 43 x 38 x 15.4cm **Pricing Unit:** unit/1 **Carton QTY:** 1 This sturdy and rugged water and airtight case is ideal for exceptionally humid or rugged conditions. This customisable case is essential for all marine and outdoor applications and is supplied with a rechargeable moisture absorbant to remove any moisture that may be present in the atmosphere after opening the case. Rated to IP67 this case can also be equipped with the Dessicant reusable Silica Gel. 



### More Information

<b>MANUFACTURER_MS</b>	Germany
<b>MODEL NUMBER_MS</b>	K61SD521
<b>ECOLABEL_MS</b>	yes