



#0 Child Disposable Guedel Airway 60mm (Black)

€2.53

In stock
SKU GAB1000

 Add to Wish List

 Add to Compare

Airway Guedel - Disposable #0

Details

Size: 60mm **Pricing Unit:** unit/1 **Carton QTY:** 50 Guedel Airways are a critical item for any advanced trauma response to keep airways clear of obstruction when providing advanced first response treatment. Guedel Airways are available in a range of sizes which must be appropriately matched to the patient to ensure a correct fit. Colours help with ease of recognition of a suitable size. Each airway is individually packaged for best hygiene practice. Guedel Airway #0 is colour black, and is 60mm in diameter.



More Information

| | |
|------------------------|----------|
| MANUFACTURER_MS | Germany |
| MODEL NUMBER_MS | K61SD521 |
| ECOLABEL_MS | yes |