



## BRAYDEN PRO CPR Manikin (Bluetooth Connectivity)

€1,754.95

In stock  
Only **395** left  
**SKU** IM16-R

 Add to Wish List

 Add to Compare

Brayden PRO Manikin (Red Blood Circulation)



### Details

**Size:** Adult **Pricing Unit:** unit/1 **Carton QTY:** 1 

The **Brayden Pro** provides healthcare professionals with **quantitative data-driven feedback**. The Brayden Pro sends live data to the [Brayden Pro App](#) via Bluetooth to provide real-time feedback of CPR performance this data can then be reviewed and analysed to improve the quality of CPR performance and provide a numeric, quantitative assessment. Up to 6 Brayden Pro's can be connected to 1 device allowing for simultaneous training and assessment. Records can be stored and easily exported via Email for easy analysis on your Spreadsheet program.

[Click Here for our Dedicated Brayden Site.](#)





### More Information

<b>MANUFACTURER_MS</b>	Germany
<b>MODEL NUMBER_MS</b>	K61SD521
<b>ECOLABEL_MS</b>	yes