



RAPIDSTOP Non-Medicated Training Haemostatic Gauze

€8.49

In stock
SKU RSH340

 Add to Wish List

 Add to Compare

RapidStop Haemostatic Gauze



Details

RapidStop Haemostatic Gauze joins the popular range of RapidStop products as a proven training device for simulation-based training. Z-folded for easy unravelling, in a vacuum sealed package, the gauze is designed for true realism. Suitable for use with the range of TraumaSim products. RapidStop Haemostatic Gauze is designed to improve student retention in



training of excessive bleeding and wound packing.



More Information

MANUFACTURER_MS	Germany
MODEL NUMBER_MS	K61SD521
ECOLABEL_MS	yes