



Wooden Non-Sterile Tongue Depressors Box/100

€4.59

In stock
SKU TDNS001

 Add to Wish List

 Add to Compare

Wooden Tongue Depressors



Details

Size: 15cm **Pricing Unit:** box/100 **Carton QTY:** 50 The AERODIAGNOSTIC range offers a selection of commonly used diagnostic supplies in the healthcare industry.

Tongue Depressor

A wooden tongue depressor is used to depress the tongue when examining the throat or mouth.



More Information

MANUFACTURER_MS	Germany
MODEL NUMBER_MS	K61SD521
ECOLABEL_MS	yes