



## QV Intensive Cream

In stock  
**SKU** QV Intensive  
Cream



Product Name	Qty
<b>QV Intensive Cream - 100g</b>	<input type="text" value="0"/>
<b>QV Intensive Cream - 500g</b>	<input type="text" value="0"/>

**€7.58**

**€17.14**

 [Add to Wish List](#)

 [Add to Compare](#)

QV Intensive Cream is a rich cream scientifically designed to moisturise and help maintain moisture in extremely dry skin. The preservative-free formulation is suitable for very sensitive skin, or skin prone to sensitivity due to preservatives.

### Details

QV Intensive Cream is a rich cream scientifically designed to moisturise and help maintain moisture in extremely dry skin. The preservative-free formulation is suitable for very sensitive skin, or skin prone to sensitivity due to preservatives.

### Benefits:

- Triple-moisturising action, including 10% glycerin, for optimal skin hydration without a greasy residue
- Clinically proven 24-hour moisturisation
- Easily applied and readily absorbed into the skin without excessive rubbing
- Helps maintain the skin's natural protective barrier and lock in moisture



- Preservative-free, further reducing the chance of irritancy
- Low-irritant formulation, fragrance, lanolin, propylene glycol and colour-free
- Non-comedogenic, so it won't block pores
- Hypoallergenic & Dermatologically tested
- Vegan Friendly - contains no animal-derived ingredients
- Cream suitable for use with sensitive and eczema-prone skin

## Directions:

Apply Intensive Cream to affected areas as often as required, particularly after washing, showering, bathing, shaving, exposure to harsh climatic conditions and at night.

## Cautions:

Avoid eyes. For external use only. Store below 30°C.

## Ingredients:

Aqua (Water), Petrolatum, Paraffinum Liquidum, Glycerin, Cyclopentasiloxane, Cyclohexasiloxane, Cetearyl Alcohol, Behentrimonium Methosulfate, Dimethicone, 1,2-Hexanediol, Caprylyl Glycol.