



WM1 Workplace First Aid Kit Wallmount - Metal Case by Trafalgar

€140.32

In stock
Only **997** left
SKU 876478

 Add to Wish List

 Add to Compare

This sturdy mounted workplace first aid get has internal shelving for quick access.



Details

Wall Mount Workplace First Aid Kit: Metal Case

Keep your workplace safe and prepared with a mounted first aid kit!

This comprehensive metal wall mounted first aid kit has your workplace covered for cuts and abrasions, sprains, fractures and bruising.

This sturdy mounted first aid kit stores all products in internal shelving for quick and easy access in times of emergency.

Dimensions:

240 x 385 x 115mm (W x H x D)



Contents

Cuts and Abrasions

- Adhesive Premium Fabric Strips x50
- Combine Dressing 10cm x 20 cm x1
- Gauze Pieces (sterile) 7.5cm x 7.5cm pack of 5 x3
- Non-Adherent Dressing 5cm x 5cm x6
- Non-Adherent Dressing 7.5cm x 10cm x3
- Non-Adherent Dressing 10cm x 10cm x1
- Wound Dressing No. 14 (medium) x1
- Wound Dressing No. 15 (large) x1
- Wound Wipes x10

Sprains, Fractures, Bruising

- Bandage - Conforming 5cm x 1.8m x3
- Bandage - Conforming 7.5cm x 1.8m x3
- Bandage - Crepe Medium Support 10cm x2m x1
- Instant Cold Pack x1
- Triangular Bandage Disposable 110cm x 155cm x2