



3M 2560 Red Dot Monitoring Electrodes (50) with Foam Tape, Sticky Gel & Press Stud



€26.00

In stock
Only **998** left
SKU 400265

 Add to Wish List
 Add to Compare



Details

Description

Clinically-proven to last up to 5 days of wear time, the 3M 2560 Red Dot Monitoring Electrodes are designed for high adhesion purposes and requirements. It is created with sticky gel to ensure excellent tracing quality for better patient information.

Made with soft foam, it feels much more comfortable than the usual electrodes to the patient. Apart from superior comfort, these electrodes also have a 45-day out of the bag performance. It is also occlusive to fluids, and easy to handle for quicker application.

Specifications

Suitable for Emergency Room (ER), Operating Room (OR), Intensive Care Unit (ICU), and Paramedics.

Electrode Type: Diaphordic, Fluid Resistant, and Monitoring

Electrode Size: 1.6" x 1.63"; 3.5cm x 4cm



Gel Type: Sticky

Skin Type: Normal

Radiolucency: Radiolucent and Radiopaque

Quantity: PKT/50

Important

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other qualified healthcare provider. Consult your healthcare provider immediately if you think that you may have a medical emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified healthcare provider before starting any new treatment or with any questions that you may have regarding a medical condition or injury.

Related Products

Check items to add to the cart or select all