



## Aluminium Finger Splint 2.5 x 46cm

€8.81

In stock  
Only **1000** left  
**SKU** 9115-04

 Add to Wish List

 Add to Compare

### Details

- The Aluminium Finger Splint is indicated to transverse fractures of the phalanges. They are curved to conform to finger in semi-flexed position. Narrow in width, this splint has an upward curve at the tip to protect the end of the finger.

This splint is not adjustable.

- Size 2.5 x 46cm

### Important

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other qualified healthcare provider. Consult your healthcare provider immediately if you think that you may have a medical emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified healthcare provider before starting any new treatment or with any questions that you may have regarding a medical condition or injury.

### Related Products

Check items to add to the cart or select all