



Basic Dressing Pack with Plastic Tray, 3 Forceps, Sterile Field, Towel & 5 Gauze

€2.31

In stock
Only **1000** left
SKU 808

 Add to Wish List

 Add to Compare

Details

Basic dressing packs are used for basic dressing and tissue procedures.

Popular for use with suture / ulcer & bed sores. Sage woven dressing packs are popular for fluid handling procedures. They contain woven swabs which generally absorb fluids at a faster rate than non-woven swabs.

- 1 x Standard Tray 145 x 130 x 18mm
- 5 x Gauze Swabs 7.5 x 7.5cm 8-Ply
- 3 x Forceps Green
- 1 x Towel 2-Ply Paper 42 x 40cm
- 1 x Field White Plastic 50 x 50cm

Important

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other qualified healthcare provider. Consult your healthcare provider immediately if you think that you may have a medical emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified healthcare provider before starting any new treatment or with any questions that you may have regarding a medical condition or injury.



Related Products

Check items to add to the cart or select all