



## I-Gel Airway

€33.93

In stock  
Only **1000** left  
**SKU** WEB-IGA

 Add to Wish List

 Add to Compare

### Details

An evolutionary supraglottic airway with non-inflatable cuff designed to mirror the patient's anatomy. The i-gel works in harmony with the patient's anatomy so that compression and displacement trauma are significantly reduced or eliminated.

### Sizes:

- #1 - Neonate
- #1.5 - Small Paediatric
- #2 - Infant
- #2.5 - Large Paediatric
- #3 - Small Adult
- #4 - Medium Adult
- #5 - Large Adult

### Important

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other qualified healthcare provider. Consult your healthcare provider immediately if you think that you may have a medical emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified healthcare provider before starting any new treatment or with any questions that you may have regarding a medical condition or injury.



---

**Related Products**

Check items to add to the cart or select all