



## Jelonet Paraffin Gauze Dressing

€19.52

In stock  
Only **1000** left  
**SKU** WEB-JPGD

 Add to Wish List

 Add to Compare



### Details

JELONET is a paraffin gauze dressing.

### Features:

- Soft paraffin base
- Sterile leno weave presentation
- Comprehensive size range

### Benefits:

- Soothes and protects the wound and allows free passage of viscous exudate.
- Jelonet is soothing and low-adherent and allows the wound to drain freely into an absorbent secondary dressing.

### Maintains shape, resists fraying

- Jelonet is a sterile paraffin tulle gras dressing made from open weave gauze. The gauze has interlocking threads which minimise fraying when the dressing is cut to shape.

### Can be combined with topical medication of choice

- Jelonet is not medicated and so is ideal for use with the topical antiseptic or antibiotic of choice

### Indications:

- Minor burns and scalds
- Donor and recipient graft sites



- Skin loss wounds, lacerations and abrasions
- Leg ulcers

**Options:**

10 x 10cm (10)

5 x 5cm (50)

**Important**

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other qualified healthcare provider. Consult your healthcare provider immediately if you think that you may have a medical emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified healthcare provider before starting any new treatment or with any questions that you may have regarding a medical condition or injury.

**Related Products**

Check items to add to the cart or select all