



## Kendrick Traction Device - KTD

€404.88

In stock  
Only **1000** left  
**SKU** 280012

 Add to Wish List  
 Add to Compare

### Details

Cost effective, compact and lightweight, the Kendrick Traction Device (KTD) is ideal for remote area operations and outdoor activities.

The traction pole length quickly adjusts for adult or paediatric application, eliminating the need for patient roll-over or unnecessary leg raising. There is no ischial bar, which reduces pressure problems.

The KTD comes complete within a very compact storage pouch with easy-to-follow instructions.

- Dimensions (Stored): 24cm l x 9cm w
- Weight: 567g

### Important

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other qualified healthcare provider. Consult your healthcare provider immediately if you think that you may have a medical emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified healthcare provider before starting any new treatment or with any questions that you may have regarding a medical condition or injury.

### Related Products

Check items to add to the cart or select all