




Murine Revital Eyes 15ml

\$25.96

In stock
Only **1000** left
SKU 409

 [Add to Wish List](#)

 [Add to Compare](#)

Details



Revital Eyes is a gently buffered solution specially formulated to closely match the natural fluid of the eye. This "natural tears formula" contains many of the major ingredients found in natural tears.

It helps relieve irritation from:

- Dust
- Air pollutants
- Sun glare
- Computer stare
- Smog
- Smoke
- Air conditioning
- Ingredients:
- Poly vinyl alcohol 14mg/ml
- Povidone 6mg/ml
- Potassium chloride
- Sodium chloride
- Sodium bicarbonate
- Sodium citrate
- Glucose
- Sodium phosphate (monobasic and dibasic) in purified water with Disodium Edetate 0.5mg/ml and Benzalkonium chloride 0.1mg/ml as preservatives.

Important

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other qualified healthcare provider. Consult your healthcare provider immediately if you think that you may have a medical emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified healthcare provider before starting any new treatment or with any questions that you may have regarding a medical condition or injury.

Related Products

Check items to add to the cart or select all