



Pacific Medical Fingertip Pulse Oximeter

\$45.87

Regular Price \$62.94

In stock
Only **80** left
SKU PACPMFPO-
TM

 Add to Wish List

 Add to Compare

This pulse oximeter is suitable for the monitoring of adult and paediatric patients in medical clinics, hospitals and at home.



Details

Stylish fingertip pulse oximeter that measures and displays SpO2, Pulse Rate and Perfusion Index

This pulse oximeter is suitable for the monitoring of adult and paediatric patients in medical clinics, hospitals and at home.

- Special splash proof and drop resistant design
- OLED display of SpO2, PR, PI, pulse bar & waveform
- Spot check and continuous pulse oximeter measurement modes
- Pulse rate analysis for spot check measurement
- Up to 12 groups SpO2 data storage
- Audible & visual alarm with low battery indication
- Automatic power on/off
- Four direction display
- Artefact removal & anti motion
- Accurate & sensitive



Pulse Oximeter Specifications:

Spo2

- Transducer: Dual-wavelength LED sensor
- Measuring range: 35-100%
- Measuring accuracy: less than or equal to 2% range from 70-100%

Pulse Rate

- Measuring range: 30-240bpm
- Measuring accuracy: plus or minus 2bpm or plus or minus 2% (whichever is the greater)

Perfusion Index

- Display range: 0-20%

Default Alarm Limit

- SpO2 low limit: 90%
- Pulse rate high limit: 120bpm
- Pulse rate low limit: 50bpm

Measuring Mode

- 30 seconds spot check mode
- Continuous mode

Dimensions, Weight & Power

- Dimensions: 56x34x30mm (L x W x H)
- Weight: 52g
- Power: 2 x AAA alkaline battery

Important

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other qualified healthcare provider. Consult your healthcare provider immediately if you think that you may have a medical emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified healthcare provider before starting any new treatment or with any questions that you may have regarding a medical condition or injury.



More Information

SPECIAL PRICE

\$45.87

Frequently bought together

You may also need



[Liberty Pocket Otoscope & Ophthalmoscope Diagnostic Set](#)

\$445.65

Add to cart