





## Red Bag with Shoulder Strap Empty 350l x 340w x 160d

**\$153.94**

In stock  
Only **1000** left  
**SKU** 350000

 Add to Wish List  
 Add to Compare



### Details

- A red lightweight paramedic's bag with a shoulder strap and fluoro stripe across the front. It has numerous compartments to store all first aid items.
  - Dimensions: 350 x 340 x 160

### Important

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other qualified healthcare provider. Consult your healthcare provider immediately if you think that you may have a medical emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified healthcare provider before starting any new treatment or with any questions that you may have regarding a medical condition or injury.

### Related Products

Check items to add to the cart or select all