



Sharps Container 4.5L Square with Flange and Seal

\$17.40

In stock
Only **999** left
SKU 9971

 Add to Wish List
 Add to Compare

Details

Sharps Container 2L Square with Flange and Seal

Type - Inter-stacking (nesting), seamless, injection molded

Capacity to fill line - 4L

Capacity to spill line - 4.75L

Lid aperture - 90mm

Final Closure (red) - 90mm

Finger guard (yellow) - 90mm

Ideal application - Bedside, trolley clip-on or wall mounting; sharps, waste, eg needles, broken glass, hypodermics, swabs

Shape - Square; with taper top to base

Dimensions - 205 d x 185 w x 185 h mm

Important

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other qualified healthcare provider. Consult your healthcare provider immediately if you think that you may have a medical emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified



healthcare provider before starting any new treatment or with any questions that you may have regarding a medical condition or injury.