



## Sqwincher 2L Liquid Concentrate

**\$52.76**

In stock  
Only **999** left  
**SKU** WEB-SQW

 Add to Wish List  
 Add to Compare

### Details

- Sqwincher® is an Electrolyte enhanced beverage for effective hydration.

Suitable for use at Water coolers and Hydration stations. Used during outdoor work, heavy machinery operations, site work, prolonged physical activity, hot weather, work with heavy PPE clothing.

- 2L bottle mixes to 20L of Sqwincher®.
- Low in sodium.
- Caffeine free.
- Flavours: Lemon Lime, Orange, Wild Berry, Tropical Cooler, Grape and Lemonade.
- Great taste for increased fluid consumption.

### Important

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other qualified healthcare provider. Consult your healthcare provider immediately if you think that you may have a medical emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified healthcare provider before starting any new treatment or with any questions that you may have regarding a medical condition or injury.

### Related Products



---

Check items to add to the cart or select all