



## Yellow Biohazard Waste Bags 560 x 990mm (200) 55L

**\$162.37**

In stock  
Only **1000** left  
**SKU** 99427

 [Add to Wish List](#)

 [Add to Compare](#)

### Details

Sentry Yellow Biohazard Waste Bags are made from a mixture of low linear density polyethylene, high density polyethylene and ethyl-butylene-polyethylene. This combination provides more strength than a bag utilising only one of these polyethylenes of the same thickness. The bag contains organic colouring and no toxic by-products are incineration.

Measurements:

565mm x 990mm

55L

50um

Quantity: 200

### Important

You must not rely on the information on this site as an alternative to medical advice from your medical physician or other qualified healthcare provider. Consult your healthcare provider immediately if you think that you may have a medical emergency or for the diagnosis and treatment of injuries. Always seek the advice of your physician or other qualified healthcare provider before starting any new treatment or with any questions that you may have regarding a medical condition or injury.

### Related Products

Check items to add to the cart or select all