



Zoll AED Plus Defibrillator (Semi/Fully Automatic)

\$3,468.17

In stock
Only **390** left
SKU 400357

Select Model

-- Please Select -- ▼

♥ Add to Wish List

■ Add to Compare

Includes:

- LCD screen showing voice prompt messages with step by step directions
- 10 x Type 123 Lithium Batteries (PN: 8000-0807-01. Five (5) year shelf-life
- CPR-D padz one piece elctrode pad with Real CPR Help (PN: 8900-0800-01). Give (5) year shelf life.
- CPR-D accessory kit (Rescue kit) including CPR barrier mask, scissors, gloves, prep razor, towel and a moist towelette.
- Soft Carry Case
- Operator's Guide and basic Life Support Poster
- Access to ZOLL Online Learning
- (7) Seven year limited warranty

This item is GST Free.

Details

• AED Plus

When it comes to helping a victim of sudden cardiac arrest, not all AEDs are created equal. Unlike other AEDs, ZOLL® AEDs support rescuers with integrated, real-time feedback on compression rate and depth. Recommended by medical professionals, the ZOLL AED Plus® defibrillator with [Real CPR Help®](#) technology guides rescuers in performing high-quality CPR and provides integrated audio and visual feedback to support you throughout the rescue.



Smart Technology

Lives are not saved by AEDs alone. People save lives — using AEDs. Real CPR Help guides you in delivering high-quality CPR. Integrated, real-time feedback provides visual and audio prompts to "Push Harder" or confirms when you are doing "Good Compressions," keeping you on track to administer potentially life-saving CPR.

In addition to step-by-step audio prompts, Real CPR Help technology includes a compression depth indicator that provides real-time visual guidance to help you deliver the best possible CPR.

Ready for the Rescue

Your AED needs to be ready to use at a moment's notice. The AED Plus defibrillator can withstand many environmental conditions and has an automatic self-test capability, so you can be sure it's ready for the rescue. And an industry-leading five-year shelf life for batteries and pads means fewer replacements, less maintenance, and a low total cost of ownership.

Trusted Partner

For more than 25 years, ZOLL has been a pioneer in resuscitation and defibrillation, designing products that push the boundaries of technology to empower first responders and unexpected heroes to help save lives.



AED Plus sets the benchmark for design, dependability, durability, and performance.

Product Features

- Real CPR Help provides real-time CPR feedback on the rate and depth of chest compressions
- Audio prompts and visual icons complement the Chain of Survival
- Easy-to-place, one-piece CPR-D-padz[®] electrodes with a 5-year shelf life
- Fully automatic feature available — unit will deliver a shock on its own when a shockable heart rhythm is detected
- Rugged design with an IP (ingress protection) rating of IP55 for dust and water
- Can withstand conditions like extreme temperature, altitude, vibration, and shock
- Off-the-shelf lithium batteries (Duracell 123a) with a 5-year battery life
- Intelligent paediatric capability: Pedi-padz[®] paediatric electrodes can be purchased separately for children up to 8 years old or weighing less than 25 Kg's
- Automatic self-testing helps ensure a state of readiness
- Soft carry case for easy travel
- 5-year warranty (7-year warranty with online registration)

"I was fumbling as I turned on the AED and applied the pads. Once I got the pads on, the AED took over. I was petrified. It was a comfort that the AED Plus walked me through the rescue ... until the emergency response team was able to take over."



Carl Summers

Blain's Farm & Fleet, Ottawa, Illinois



More Information

BRAND

ZOLL

Frequently bought together

You may also need



[Zoll Adult CPR-D Pads for AED Plus & Pro - Depth & Rate Sensors](#)

\$230.43

Add to cart



[Zoll AED Plus Trainer II Complete](#)

\$621.94

[Add to cart](#)



[Zoll Soft Carry Case for AED Plus & AED Pro](#)

\$245.36

[Add to cart](#)



[Zoll AED Plus Alarmed Wall Mounting Box & Surface Unit](#)

\$662.48

[Add to cart](#)