



QV Bar Twin Pack - 2x 100g


\$18.66

In stock
Only **1000** left
SKU 10090

 Add to Wish List

 Add to Compare

The QV Bar that looks like soap. It feels like soap. But it's not soap. It's better.



Details

The QV Bar that looks like soap. It feels like soap. But it's not soap. It's better.. QV Bar is enriched with moisturisers, and will leave your skin feeling soft and smooth.

Directions:

Wet skin, apply QV Bar to desired areas. Rinse pat skin gently dry. Do not use soap. Keep bar dry for prolonged usage. For persistent dry skin use in conjunction with QV Bath Oil and QV Cream. Tip - Always allow your bar to drain well between uses to get the most out of your product.

Cautions:

Avoid eyes, if lather enters eyes, wash out immediately with water. For external use only. Store below 30°C.



Ingredients:

Disodium Laureth Sulfosuccinate, Cetearyl Alcohol, Sodium Lauryl Sulfate, Disodium Cocamido MEA-Sulfosuccinate, PEG-150, ZEA Mays (Corn) Startch, Stearic Acid, Cocamidopropyl Betaine, Glycerin, Petrolatum, Stearyl Heptanoate, Dimethicone, Paraffin, Lactic Acid, Methylparaben.