



QV Intensive Moisturiser - 100g

\$9.24

In stock
Only **500** left
SKU 10180

 Add to Wish List

 Add to Compare

A highly emollient water-free soothing ointment for extremely dry, cracked sensitive skin. Helps to rehydrate and seal moisture in the the skin.



Details

QV Intensive Moisturiser is a highly emollient balm that has been formulated to care for extremely dry or cracked skin. The moisturising balm helps to seal moisture into the skin, and the water-free formulation won't sting on cracked skin.

Benefits:

- Water-free formulation won't sting when applied to inflamed or cracked skin.
- Suitable for dry, sensitive, itchy skin and eczema.
- Low-irritant formulation, free from fragrance, colour and lanolin to suit sensitive skin types.
- Helps to maintain skin's protective properties.
- Vegan Friendly - contains no animal derived ingredients.
- Body moisturising cream suitable for use with sensitive and eczema prone skin.



Directions:

Apply body moisturiser liberally to the face and body and massage in gently. Use after bathing and showering.

Cautions:

Avoid eyes. For external use only. Store below 30°C.

Ingredients:

Light Liquid Paraffin (Paraffinum Liquidum), Petrolatum, Isopropyl Myristate, Polyethylene, Cetearyl Alcohol, Silica, Dimethicone