



AEROFORM Finger Elastic Tubular Bandage Size 01 x 20M

\$100.63

In stock
SKU AB0120

 [Add to Wish List](#)

 [Add to Compare](#)

Size 01 Finger Tubular Bandage 20M

Details

Size: Size 01 x 20M **Pricing Unit:** unit/1 **Carton QTY:** 6 AEROFORM Tubular Bandages are an elasticised support bandage suited to treatment of sports injuries, dislocations and sprains.

Self Supporting

AEROFORM Tubular Bandages are held in place without tying or tapes, resulting in a comfortable long wearing bandage.

Provides Even Support

AEROFORM Tubular Bandages provide even support across damaged joints and muscles due to its elastic self supporting properties resulting in a comfortable support bandage.



More Information

MANUFACTURER_MS	Germany
MODEL NUMBER_MS	K61SD521
ECOLABEL_MS	yes