



AEROCREPE Elastic Crepe Bandage 15cm x 4M Wrap/12

\$80.26

In stock
SKU ACE150

 Add to Wish List

 Add to Compare

Elastic Crepe Bandage 15cm


Details

Size: 15cm x 4M **Pricing Unit:** wrap/12 **Carton QTY:** 20 The AEROCREPE Bandage is a high quality cotton crepe bandage, recommended for low support requirements and can be used for the retention of dressings.

Elastic Crepe

AEROCREPE elastic cotton crepe bandages are of an elastic conforming nature which is durable and suitable for immobilisation and compression techniques.

Easy to Open Wrapper

AEROCREPE elastic cotton crepe bandages have a easy open tab built into the wrapper of the bandage. 



More Information

| | |
|------------------------|----------|
| MANUFACTURER_MS | Germany |
| MODEL NUMBER_MS | K61SD521 |
| ECOLABEL_MS | yes |