



## AEROCREPE Elastic Crepe Bandage 7.5cm x 4M Wrap/12

**\$43.78**

In stock  
**SKU** ACE75

 Add to Wish List

 Add to Compare

Elastic Crepe Bandage 7.5cm


### Details

**Size:** 7.5cm x 4M **Pricing Unit:** wrap/12 **Carton QTY:** 30 The AEROCREPE Bandage is a high quality cotton crepe bandage, recommended for low support requirements and can be used for the retention of dressings.

### Elastic Crepe

AEROCREPE elastic cotton crepe bandages are of an elastic conforming nature which is durable and suitable for immobilisation and compression techniques.

### Easy to Open Wrapper

AEROCREPE elastic cotton crepe bandages have a easy open tab built into the wrapper of the bandage. 



## More Information

<b>MANUFACTURER_MS</b>	Germany
<b>MODEL NUMBER_MS</b>	K61SD521
<b>ECOLABEL_MS</b>	yes