



## AEROCREPE Light Cotton Crepe Bandage 15cm x 4M Wrap/12

**\$88.79**

In stock  
**SKU** ACL150

 Add to Wish List

 Add to Compare

Light Crepe Bandage 15cm




### Details

**Size:** 15cm x 4M **Pricing Unit:** wrap/12 **Carton QTY:** 20 The AEROCREPE Bandage is a high quality cotton crepe bandage, recommended for low support requirements and can be used for the retention of dressings.

### Light Crepe

AEROCREPE light cotton crepe bandages are of a standard conforming nature which is durable and suitable for dressing retention and light support requirements.

### Easy to Open Wrapper

AEROCREPE light cotton crepe bandages have a easy open tab built into the wrapper of the bandage. 



### More Information

<b>MANUFACTURER_MS</b>	Germany
<b>MODEL NUMBER_MS</b>	K61SD521
<b>ECOLABEL_MS</b>	yes